



## **OPMI PENTERO 900**

The Next Generation



We make it visible.





# OPMI PENTERO 900

## The Next Generation

### Evolution in Excellence

OPMI® PENTERO® 900 represents the next generation visualization system. Building on the groundbreaking innovations introduced in 2004, it combines unique design concepts and new functionalities in a proven, fully-integrated platform. Key functions have been enhanced and new visualization methods integrated, raising OPMI PENTERO 900 to a new level of performance. OPMI PENTERO 900 continues the evolutionary process, turning technological progress into medical innovations effectively advancing what is possible in modern microsurgery.

Experience the new OPMI PENTERO 900.

### ■ Brilliant Visualization


Experience optical immersion with state-of-the-art apochromatic optics, razor-sharp video images in full HD quality and innovative fluorescence methods.

### ■ Superior Performance

Smooth, intuitive system handling and superior functionality ensure efficient surgery and fast system set-up.

### ■ Beyond Visualization

OPMI PENTERO 900 interacts with current and emerging workplace technologies and workflow based solutions to create a better OR experience.



// INNOVATION  
MADE BY CARL ZEISS



# // BRILLIANT VISUALIZATION

MADE BY CARL ZEISS

## Experience Optical Immersion

OPMI PENTERO 900 delivers state-of-the-art apochromatic optics providing crystal-clear images, sharp details and natural color rendition. Whether viewing through the eyepieces or on its flexible and integrated HD touchscreen display, OPMI PENTERO 900 elevates visualization of the surgical field to the next level. The entire HD video chain – camera, recorder, editor and monitor – is fully integrated into the system without the need for external components, exposed cabling or the use of multiple control interfaces. The HD video system can be configured and controlled via the central HD touchscreen, handgrip or foot switch for maximum flexibility and surgical performance. OPMI PENTERO 900 offers a unique HD experience with visual brilliance for live demonstrations, teaching presentations and patient documentation.



**Uniquely Designed Apochromatic Optics**  
The distinctive design concept with ZEISS apochromatic optics throughout the entire optical pathway allows the system to deliver unmatched optical clarity, detail resolution and color reproduction, both through the eyepieces and the video image.



**Fully Integrated HD Video**  
The fully integrated HD video camera, recorder and editor enable surgeons to capture razor-sharp images for teaching, documentation and presentation purposes. All video functions can be centrally controlled from the intuitive graphical user interface.




**Brilliant High Definition Video Display**  
The large touchscreen display delivers impressively, crisp images in HD quality. The extendable suspension arm allows the display to be rotated, tilted or moved into different viewing positions.




**High-quality Data Injection**  
The advanced MultiVision™ display significantly enhances image quality and enables efficient data injection during procedures. Higher resolution, enhanced contrast and better color rendition ensure outstanding quality leading to improved outcomes.

## Continuation of Scientific Progress


Carl Zeiss has partnered with visionary thought leaders in the development of innovative intraoperative fluorescence technologies. The first fluorescence module, INFRARED 800™, successfully established intraoperative fluorescence as a key visualization technology. FLOW® 800 is a unique tool enabling the visual analysis of blood flow dynamics, further establishing Carl Zeiss as a leader in intraoperative fluorescence. With BLUE 400 and now, with the introduction of YELLOW 560™, we further broaden the capability for fluorescencebased research applications. With its intuitive workflow, automated functions and unmatched performance, the PENTERO 900 platform supports fluorescence-based surgeries like no other system.




**INFRARED 800<sup>1</sup>**  
Intraoperative visual assessment of blood flow and vessel patency during AVM, bypass and aneurysm surgery. INFRARED 800 is indicated for use in neurosurgery, plastic and reconstructive procedures and coronary artery bypass graft surgery.



**FLOW 800<sup>1</sup>**  
Unique fluorescence application enabling visual analysis of vascular blood flow dynamics. It compiles INFRARED 800 video sequences into visual maps, diagrams or side-by-side images, enabling an in-depth interpretation of fluorescence videos.



**BLUE 400<sup>2</sup>**  
Capable of supporting fluorescence-based surgery for research applications by providing visualization in the 620-710 nm range. BLUE 400 is a completely integrated module with no external components or cables. Optionally available in HD quality.



**YELLOW 560<sup>3</sup>**  
Visualizes in the 540 to 690 nm wavelength range for supporting research applications. It is the first ZEISS fluorescence module capable of highlighting fluorescent structures while viewing surrounding structures in their natural color.

<sup>1</sup> For the complete 510(k) summary for INFRARED 800 with FLOW 800 option, visit [http://www.accessdata.fda.gov/cdrh\\_docs/pdf10/k100468.pdf](http://www.accessdata.fda.gov/cdrh_docs/pdf10/k100468.pdf).

<sup>2</sup> Image obtained with investigational new drug.

<sup>3</sup> Image obtained with drug under investigation for new intended use.



# // SUPERIOR PERFORMANCE MADE BY CARL ZEISS

## Maximum Efficiency

OPMI PENTERO 900 is a unique surgical visualization platform specifically conceived and designed for even the most demanding microsurgical applications. Extended system ergonomics and functions provide increased convenience, streamlining the surgical workflow. All relevant functions are combined into a cohesive system that can be controlled from the intuitive touchscreen user interface. Smooth system handling and superior performance are delivered through proven functions like mouth switch control, a unique depth of field, AutoFocus™ and many more. The Foldable Tube f170/f260 and wireless foot control panel enhance ergonomics, improving surgeon comfort and performance.

## OR Turn-Around Simplified

Designed as much for the OR staff as for the surgeon, OPMI PENTERO 900 incorporates workflow-conducive features that reduce prep time for the nursing staff before each surgical case. AutoBalance™ quickly balances the microscope at the touch of a button, AutoDrape® facilitates a quick and easy draping process, and the unique FlexiTrak™ enables the OR staff to easily maneuver the system in the clinical environment. Additionally, the intuitive user interface allows for easy access to all microscope functions through a large, HD touchscreen monitor including patient data, pre-configured surgeon settings and video recording. More than any other surgical microscope system, OPMI PENTERO 900 streamlines the surgical workflow and maximizes OR efficiency.




**Dynamic Foldable Tube**  
The highly flexible Foldable Tube f170/f260 offers increased positioning capability, magnification and user comfort. With a quick dial turn, the integrated PROMAG™ functionality provides an additional 50 percent magnification gain.



**Unparalleled Depth of Field**  
Depending on preference or application, the integrated, electronically controlled double-iris diaphragm enables surgeons to choose between maximum light and resolution or depth of field.



**Efficient Light Settings**  
The patented, two-channel illumination design reduces shadowing in deep cavities. The Automatic Iris Control™ limits illumination to the field of interest and the new Focus Light Link™ automatically limits brightness, both in the aim of preventing inadvertent light exposure.



**Wireless Foot Control**  
The wireless foot control, designed to manage multiple microscope functions, can be freely positioned for rapid set-up and user convenience during surgery. The intelligent power management function ensures superior long-lasting operability.



**AutoDrape**  
AutoDrape allows for easy and fast system draping by automatically evacuating air from the sterile drape. The drape clings tightly to the system in a matter of seconds without obstructing the microscope's mobility.



**AutoBalance**  
Automatic balancing of the microscope with the press of a single button for fast system set-up. AutoBalance can even be performed while the system is draped ensuring optimum system mobility. A unique fine balance function also allows precise system control via the mouth switch.



**VisionGuard Drape Technology**  
ZEISS drapes are manufactured with VisionGuard®, a unique, optical lens that works with the microscope's objective as a single optical unit for unparalleled optical clarity. The lens can be replaced with a new sterile lens during surgery without compromising sterility.



**Rapid Remote Diagnosis**  
Fast internet access via VPN enables remote diagnostic capability, resulting in better service support and improved system uptime. System log files can be accessed online by ZEISS service specialists for rapid analysis and support.







OPMI PENTERO 900

Standard and Optional Features



Standard Configuration

Apochromatic Optics	Motorized focus; Varioskop® with working distance 200 – 500 mm
	Motorized zoom; 1:6 zoom ratio
	10x magnetic widefield eyepieces with integrated eyecups
	Autofocus with 2 visible laser dots, automatic mode with magnetic brakes
Illumination	Superlux® 330 light source with 2 x 300 W xenon
	Automatic Iris Control for adjusting the illumination to the field of view
	Individual light threshold setting
	Focus Light Link: working distance controlled light intensity
	Display of remaining lamp life time
System Operation	Multifunctional programmable handgrips
	Magnetic clutches for all system axes
	Central user interface
	XY robotic movement in 3 axes (variable speed)
System Setup	AutoBalance
	AutoDrape – air evacuation system
	Mouthswitch fine balance
Video	22" HD video touchscreen on extendable arm
	Integrated 3-CMOS SD video camera
	Integrated video still image capturing on HDD and USB-media
Connectivity / Data Management	Video-in for external SD video sources
	Navigation interface
	Interface for micromanipulator
	Remote diagnosis via internet/VPN

Options

Video	Integrated 3-CMOS HD video camera
	Integrated SD or HD video recording and editing
	Adaptation of consumer (SLR) photo/video camera
Intraoperative Fluorescence	BLUE 400
	INFRARED 800
	INFRARED 800 with FLOW 800 YELLOW 560
Connectivity / Data Management	DICOM module for patient data transfer from/to PACS
	Accessories
Accessories	12.5x magnetic widefield eyepieces with integrated eyecups
	Straight tube, focal length f = 170mm
	Stereo co-observation tube
	Foldable Tube f170/f260, including the PROMAG function for additional 50% magnification and integrated rotate function
	Wired foot control panel
	Wireless foot control panel

Image Courtesy of:

BrainLAB AG, Feldkirchen, Germany (p. 4, 8)

Barrow Neurological Institute, Phoenix, Arizona, USA (p. 1, 4, 5, 8)

Michael Buchfelder, MD, PhD, Neurosurgery Department, Erlangen University Hospital, Erlangen, Germany (p. 5)

Walter Stummer, MD, PhD, Department of Neurosurgery, University Hospital Münster, Münster, Germany (p. 5, 8)

Yasuo Murai, MD, PhD, Department of Neurosurgery, Nippon Medical School, Tokyo, Japan (p. 8)

Yasushi Takagi, MD, PhD, Department of Neurosurgery, Kyoto University Graduate School of Medicine, Kyoto, Japan (p. 5)



#### Your local contact:

##### Argentina

Carl Zeiss Argentina S.A.  
Calle Nahuel Huapi 4015 / 25  
C1430 BCO Buenos Aires  
Argentina  
Phone: +54 11 45 45 66 61  
bruzzi@zeiss.com.ar

##### Australia

Carl Zeiss Pty. Ltd.  
Unit 13, 2 Eden Park Drive  
North Ryde, New South Wales 2113  
Australia  
Phone: +61 2 9020 1333  
med@zeiss.com.au

##### Austria

Carl Zeiss GmbH  
Laxenburger Str. 2  
1100 Vienna  
Austria  
Phone: +43 1 79 51 80  
austria@zeiss.org

##### Belgium

Carl Zeiss NV-SA  
Ikaroslaan 49  
1930 Zaventem  
Belgium  
Phone: +32 2 719 39 11  
info@zeiss.be

##### Brazil

Carl Zeiss do Brasil Ltda.  
Av. Nações Unidas, 21711  
CEP04795-100 São Paulo  
Brazil  
Phone: +55 11 5693 5521  
medbrasil@zeiss.org

##### Canada

Carl Zeiss Canada Ltd.  
45 Valleybrook Drive  
Toronto, ON M3B 2S6  
Canada  
Phone: +1 800 387 8037  
micro@zeiss.com

##### China

Carl Zeiss Shanghai Co. Ltd.  
1/f., Ke Yuan Building  
11 Yi Rin Nan Road  
Waigaoqiao Free Trade Zone  
2005 Yang Gao Bei Road  
Shanghai 200131  
China  
Phone: +86 21 5048 17 17  
sro@zeiss.com.cn

##### Czech Republic

Carl Zeiss spol. s.r.o.  
Radlická 14/3201  
150 00 Prague 5  
Czech Republic  
Phone: +420 233 101 221  
zeiss@zeiss.cz

##### France

Carl Zeiss Meditec France SAS  
60, route de Sartrouville  
78230 Le Pecq  
France  
Phone: +33 1 34 80 21 00  
med@zeiss.fr

##### Germany

Carl Zeiss Meditec VG mbH  
Carl-Zeiss-Strasse 22  
73447 Oberkochen  
Germany  
Phone: +49 7364 20 6000  
vertrieb@meditec.zeiss.com  
Surgical Ophthalmology:  
Phone: +49 800 470 50 30  
iol.order@meditec.zeiss.com

##### Hong Kong

Carl Zeiss Far East Co. Ltd.  
Units 11-12, 25/F  
Tower 2, Ever Gain Plaza  
No. 88 Container Port Road  
Kwai Chung  
Hong Kong  
Phone: +852 2332 0402  
czfe@zeiss.com.hk

##### India

Carl Zeiss India Pvt. Ltd.  
22, Kensington Road  
Ulsoor  
Bangalore 560 008  
India  
Phone: +91 80 2557 88 88  
info@zeiss.co.in

##### Italy

Carl Zeiss S.p.A.  
Viale delle Industrie 20  
20020 Arese (Milan)  
Italy  
Phone: +39 02 93773 1  
post@zeiss.it

##### Japan

Carl Zeiss Meditec Japan Co. Ltd.  
Shinjuku Ku  
Tokyo 160-0003  
22 Honchio-Cho  
Japan  
Ophthalmic instruments:  
Phone: +81 3 33 55 0331  
medsales@zeiss.co.jp  
Surgical instruments:  
Phone: +81 3 33 55 0341  
cmskoho@zeiss.co.jp

##### Malaysia

Carl Zeiss Sdn Bhd.  
Lot2, Jalan 243/51 A  
46100 Petaling Jaya  
Selangor Darul Ehsan  
Malaysia  
Phone: +60 3 7877 50 58  
malaysia@zeiss.com.sg

##### Mexico

Carl Zeiss de México S.A. de C.V.  
Avenida Miguel Angel de Quevedo  
496  
04010 Mexico City  
Mexico  
Phone: +52 55 59 99 0200  
cz-mexico@zeiss.org

##### Netherlands

Carl Zeiss B.V.  
Trapezium 300  
Postbus 310  
3364 DL Slidrecht  
Netherlands  
Phone: +31 184 43 34 00  
info@zeiss.nl

##### New Zealand

Carl Zeiss (N.Z.) Ltd.  
15B Paramount Drive  
P.O. Box 121 - 1001  
Henderson, Auckland 0650  
New Zealand  
Phone: +64 9 838 5626  
med@zeiss.com.au

##### Poland

Carl Zeiss sp. z o.o.  
ul. Lopuszanska 32  
02-220 Warsaw  
Poland  
Phone: +48 22 858 2343  
medycyna@zeiss.pl

##### Singapore

Carl Zeiss Ptd. Ltd.  
50 Kaki Bukit Place  
Singapore 415926  
Singapore  
Phone: +65 6741 9600  
info@zeiss.com.sg

##### South Africa

Carl Zeiss (Pty.) Ltd.  
363 Oak Avenue  
Ferndale  
Randburg 2194  
South Africa  
Phone: +27 11 886 9510  
info@zeiss.co.za

##### South Korea

Carl Zeiss Co. Ltd.  
Seoul 121-828  
Mapo-gu  
141-1, Sangsu-dong  
2F, BR Elitel Bldg.  
South Korea  
Phone: +82 2 3140 2600  
korea@zeiss.co.kr

##### Spain

Carl Zeiss Meditec Iberia S.A.  
Ronda de Poniente, 15  
Tres Cantos  
28760 Madrid  
Spain  
Phone: +34 91 203 37 00  
info@zeiss.es

##### Sweden

Carl Zeiss AB  
Tegeluddsvaegen 76  
10254 Stockholm  
Sweden  
Phone: +46 84 59 25 00  
info@zeiss.se

##### Switzerland

Carl Zeiss AG  
Feldbachstrasse 81  
8714 Feldbach  
Switzerland  
Phone: +41 55 254 7534  
med@zeiss.ch

##### Thailand

Carl Zeiss Thailand  
Floor 8, Thosapol Land Building 2  
230 Ratchadapisek Road  
Huaykwang, Bangkok 10310  
Thailand  
Phone: +66 2 2 74 06 43  
thailand@zeiss.com.sg

##### United Kingdom

Carl Zeiss Ltd.  
15-20 Woodfield Road  
Welwyn Garden City  
Hertfordshire, AL7 1JQ  
United Kingdom  
Phone: +44 1707 871200  
info@zeiss.co.uk

##### United States of America

Carl Zeiss Meditec, Inc.  
5160 Hacienda Drive  
Dublin, CA 94568  
USA  
Phone: +1 925 557 4100  
info@meditec.zeiss.com



0297



#### Carl Zeiss Meditec AG

Goeschewitzer Strasse 51–52  
07745 Jena  
Germany  
www.meditec.zeiss.com/pentero900  
www.meditec.zeiss.com/contacts

#### Carl Zeiss Meditec Inc.

5160 Hacienda Drive  
Dublin, CA 94568  
USA  
info@meditec.zeiss.com  
www.meditec.zeiss.com

US\_30\_010\_900111

The contents of the brochure may differ from the current status of approval of the product in your country. Please contact our regional representative for more information.  
Subject to change in design and scope of delivery and as a result of ongoing technical development. OPMI, FLOW, VisionGuard, Varioskop and AutoDrape are registered trademarks of Carl Zeiss. Multivision, INFARED 800, BLUE 400, YELLOW 560, PROMAG, AutoBalance, FlexTrak, Automatic Iris Control and Focus Light Link are trademarks of Carl Zeiss. Contraves® is a registered trademark of the Contraves AG.

© 2011 by Carl Zeiss Meditec AG. All copyrights reserved. Printed in USA. 1011 2.5M